





# pcDNA<sup>™</sup>3.1/His A, B, and C

Catalog number V385-20

Revision date 2 March 2012 Publication Part number 28-0111

MAN0000627

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.



## Contents

Kit Contents and Storage	iv
Methods	1
Cloning into pcDNA <sup>™</sup> 3.1/His A, B, and C	1
Transformation and Transfection	
Appendix	9
Map of pcDNA <sup>™</sup> 3.1/His A, B, and C Vectors	
Features of pcDNA <sup>™</sup> 3.1/His A, B, and C Vectors	
Map of pcDNA <sup>™</sup> 3.1/His/ <i>lac</i> Z	11
Accessory Products	12
Technical Support	13
Purchaser Notification	14
References	15

## Kit Contents and Storage

Shipping and Storage	pcDNA <sup>™</sup> 3.1/His vectors are shipped on wet ice. Upon receipt, store vectors at –20°C.
Kit Contents	<ul> <li>10 μg each of pcDNA<sup>™</sup>3.1/His A, B, and C are supplied at 0.5 μg/μL in</li> <li>10 mM Tris-HCl, 1 mM EDTA, pH 8.0 in a total volume of 20 μL.</li> <li>10 μg of pcDNA<sup>™</sup>3.1/His/<i>lacZ</i> is supplied at 0.5 μg/μL in 10 mM Tris-HCl,</li> <li>1 mM EDTA, pH 8.0 in a total volume of 20 μL.</li> </ul>
	For research use only. Not intended for any human or animal therapeutic or diagnostic use.

### Methods

# Cloning into pcDNA<sup>™</sup>3.1/His A, B, and C

Description of the System	pcDNA <sup>TM</sup> 3.1/His A, B, and C are 5.5 kb vectors derived from pcDNA <sup>TM</sup> 3.1 and designed for high-level expression and purification of recombinant proteins in mammalian hosts. The vectors are supplied in three reading frames to facilitate in frame cloning with a polyhistidine metal-binding tag. The human cytomegalovirus (CMV) immediate-early promoter provides high-level expression in a wide range of mammalian cells. In addition, the vector replicates episomally in cell lines that are latently infected with SV40 or that express the SV40 large T antigen (e.g. COS7). High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The control plasmid, pcDNA <sup>TM</sup> 3.1/His/ <i>lacZ</i> , is the pcDNA <sup>TM</sup> 3.1/His B vector with a 3.2 kb fragment containing the $\beta$ -galactosidase gene cloned in frame with the N-terminal peptide. pcDNA <sup>TM</sup> 3.1/His/ <i>lacZ</i> is included for use as a positive control for transfection, expression, and purification in the cell line of choice.
General Molecular Biology Techniques	For help with DNA ligations, <i>E. coli</i> transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, see <i>Molecular Cloning: A Laboratory Manual</i> (Sambrook <i>et al.</i> , 1989) or <i>Current Protocols in Molecular Biology</i> (Ausubel <i>et al.</i> , 1994).
Maintaining pcDNA <sup>™</sup> 3.1/His	Many <i>E. coli</i> strains are suitable for the growth of this vector. To propagate and maintain pcDNA <sup>M</sup> 3.1/His A,B, and C, use the supplied stock solution to transform a <i>rec</i> A (recombination deficient), <i>end</i> A (endonuclease A deficient) <i>E. coli</i> strain like TOP10, TOP10F', DH5 $\alpha^{M}$ -T1 <sup>R</sup> , DH10B <sup>M</sup> , or equivalent (see page 12 for ordering information). Select the transformants on LB plates containing 50–100 µg/mL ampicillin. For long-term storage, prepare a glycerol stock of your plasmid-containing <i>E. coli</i> strain.

#### Kozak Sequence for Mammalian Expression

If you are recombining your entry clone with a destination vector for mammalian expression, your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position –3 and the G at position +4 (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined.

#### (G/A)NN<u>ATG</u>G



The pcDNA<sup>™</sup>3.1/His vectors are fusion vectors. To ensure proper expression of your recombinant protein, you must clone your gene in frame with the ATG at base pairs 920–922. This will create a fusion with the N-terminal polyhistidine tag, Xpress<sup>™</sup> epitope, and the enterokinase cleavage site. The vector is supplied in three reading frames to facilitate cloning. See pages 3–5 to develop a cloning strategy.

If you wish to clone as close as possible to the enterokinase cleavage site, follow the guidelines below:

- Digest pcDNA<sup>™</sup>3.1/His A, B, or C with Kpn I
- Create blunt ends with T4 DNA polymerase and dNTPs
- Clone your blunt-ended insert in frame with the lysine codon (AAG) of the enterokinase recognition sequence.

Following enterokinase cleavage, no vector-encoded amino acid residues will be present in your protein.

**Multiple Cloning Site of pcDNA<sup>TM</sup>3.1/His A** Below is the multiple cloning site for pcDNA<sup>TM</sup>3.1/His A. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotide indicates the variable region. Note that there is a stop codon after the *Xba* I site and that the *Asp718 I* and *Kpn I* sites are in the same reading frame for all three vectors. The multiple cloning site has been confirmed by sequencing and functional testing. The sequence is available for downloading from <u>www.lifetechnologies.com</u> or by contacting Technical Support (see page 13).

	T7 promoter/priming site								
839	CACTGCTTAC TGGCTTATCG AAATTAATAC GACTCACTAT AGGGAGACCC AAGCTGGCTA								
	Hind III Polyhistidine (6xHis) region								
899	GCGTTTAAAC TTAAGCTTAC C ATG GGG GGT TCT CAT CAT CAT CAT CAT CAT Met Gly Gly Ser His His His His His His								
950	GGT ATG GCT AGC ATG ACT GGT GGA CAG CAA ATG GGT CGG GAT CTG TAC								
	Gly Met Ala Ser Met Thr Gly Gly Gln Gln Met Gly Arg Asp Leu Tyr								
	Xpress <sup>™</sup> Epitope Asp718 I Kpn I BamH I BstX I <sup>*</sup> EcoR I								
998	GAC GAT GAC GAT AAG GTA CCT AGG ATC CAG TGT GGT GGA ATT CTG CAG Asp Asp Asp Asp Lys Val Pro Arg Ile Gln Cys Gly Gly Ile Leu Gln								
	Enterokinase recognition sequence EK cleavage site								
1040	EcoR V BstX I* Not I Xho I Xba I Apa I								
1046	ATA TCC AGC ACA GTG GCG GCC GCT CGA GTC TAG AGGGCCCGTT TAAACCCGCT Ile Ser Ser Thr Val Ala Ala Ala Arg Val ***								
	BGH reverse priming site								
1099	GATCAGCCTC GACTGTGCCT TCTAGTTGCC AGCCATCTGT TGTTTGCCCC TCCCCCGTGC								
	BGH poly (A) site								
1159	CTTCCTTGAC CCTGGAAGGT GCCACTCCCA CTGTCCTTTC CTAATAAAAT GAGGAAATTG								
*Note th	at there are two BstX I sites in the polylinker.								

Site of	<b>Balow</b> is the multiple cloning site for pcDNA <sup>™</sup> 3.1/His B. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. The multiple cloning site has been confirmed by sequencing and functional testing. Note that the <i>Asp718 I</i> and <i>Kpn I</i> sites are in the same reading frame for all three vectors. The sequence is available for downloading from www.lifetechnologies.com or by contacting Technical Support (see page 13). T7 promoter priming site									e variable d e reading from					
839	ርልርጥ	30 TT A	TGGC	ጥጥ አጥ ፖ	ים אי									ممددر	IGGCTA
055	CACIC	JUIIAC	Hind I					S GA	JICA					) region	
899	GCGTT	ΓΤΑΑΑΟ			AC C						CAT	CAT	CAT		CAT
950	GGT A Gly N	ATG GC Met Al	T AGC a Ser	ATG Met	ACT Thr	GGT Gly	GGA Gly	CAG Gln	CAA Gln	ATG Met	GGT Gly	CGG Arg	GAT Asp	CTG Leu	TAC Tyr
		Xpress™	Epitope		Asp718	SI Kpn	I Bsu	36 I B	amH I		BstX	I* <i>Eco</i>	RI		
998	Asp A	Asp As	C GAT p Asp	AAG Lys	GTA Val	CCT Pro	AAG Lys	GAT	CCA	GTG Val	TGG Trp	TGG Trp	AAT Asn	TCT Ser	GCA Ala
			gnition se					_		_				_	
1046	EcoR		.G CAC	BstX I*	I		Xho				CAC	ccc			አአአ
1040			n His												
				BC	H reve	erse pri	ning si	te							
1094			T CAG p Gln												
1142			T CCC						CCC	IGGA	AGG !	rgccz	ACTC	CC	
			BGH	poly (A	() site										
1189	ACTG	ГССТТІ	CCTA	ATAA	AA' TO	GAGG	AAAT	Г							
45.7 / A	1														

\*Note that there are two BstX I sites in the polylinker.

Multiple Cloning
 Site of
 pcDNA<sup>™</sup>3.1/His C
 Below is the multiple cloning site for pcDNA<sup>™</sup>3.1/His C. Restriction sites are labeled to indicate the cleavage site. The multiple cloning site has been confirmed by sequencing and functional testing. Note that the *Asp718 I* and *Kpn I* sites are in the same reading frame for all three vectors. The sequence is available for downloading from www.lifetechnologies.com or by contacting Technical Support (see page 13).

	T7 promoter priming site											
839	CACTGCTT	AC TGGCTTAI	CG AAA		C GAG	TCAC	CTAT	AGG	GAGAC	ccc /	AAGCI	IGGCTA
		Hind III						Polyhi	stidine	(6xHis)	) region	
899	GCGTTTAA	AC TTAAGCTI		rg GGG et Gly			-		CAT His	-		·
950		GCT AGC ATC Ala Ser Met										
	Xpress	™ Epitope	Asp718 I	Kpn I Bai	nH I		Bst.	X I* 1	EcoR I			
998	GAC GAT Asp Asp	GAC GAT AAG Asp Asp Lys		ro Gly								
	Enterokinase r	ecognition sequenc	EK clea	vage site								
	EcoR V	BstX I*	Not I	Xho I		Xba I			Apa I			
1046		GCA CAG TGG Ala Gln Trp								GTT Val	TAA ***	
		BG	H reverse pi	riming site								
1091	ACCCGCTG	AT CAGCCTCO	AC TGT	GCCTTC	r AG1	TGCC	CAGC	CATO	CTGTI	GT :	TTGCC	CCTCC
										BGH p	oly (A)	site
1151	CCCGTGCC	TT CCTTGACC	CT GGA	AGGTGC	C ACI	CCCA	ACTG	TCCI	TTCC	CTA A	ATAA	ATGAG
1211	GAAATTGC	АТ										

\*Note that there are two BstX I sites in the polylinker.

### **Transformation and Transfection**

<i>E. coli</i> Transformation	Transform your ligation mixtures into a competent <i>recA</i> , <i>endA E</i> . <i>coli</i> strain (e.g., TOP10, TOP10F', DH5 $\alpha^{\text{T}}$ -T1 <sup>R</sup> , DH10B <sup><math>\text{M}</math></sup> , see page 12) and select on LB plates containing 50–100 µg/mL ampicillin. Select 10–20 clones and analyze for the presence and orientation of your insert.					
Applying Selective Pressure	We recommend taking some (if not all) or your clone from being "overrun" by back					
	• Use carbenicillin instead of ampicillin ampicillin, and allows for a longer per preserving your clones longer.					
	• Increase the antibiotic concentration clones will not be overwhelmed by β					
	• <b>Periodically refresh plate media.</b> If y beginning to fail, spin them down, re wells with fresh LB media plus glyce	move the old media, and replenish the				
	<b>Streak clones on selective (preferably ca</b> 12 hours, isolate colonies for downstream clones from potential background contan	n usage. This will isolate your desired				
	We recommend that you sequence your of Reverse primers (see page 12 for ordering is fused in frame with the N-terminal His the diagrams on pages 3–5 for the sequer	g information) to confirm that your gene s tag and the enterokinase site. Refer to				
Ū.	Primer	Sequence				
	T7 Promoter	5'-TAATACGACTCACTATAGGG-3'				
	BGH Reverse	5'-TAGAAGGCACAGTCGAGG-3'				
	For your convenience, we offer a custom visit <u>www.lifetechnologies.com</u> or call Te	<b>1</b>				
Plasmid Preparation	Plasmid DNA for transfection into eukary from phenol and sodium chloride. Conta interfere with lipid complexing, decreasing recommend isolating DNA using the Pure PureLink <sup>®</sup> HiPure Midiprep Kit (see page gradient centrifugation.	minants will kill the cells and salt will ng transfection efficiency. We eLink® HiPure Miniprep Kit or the				

# Transformation and Transfection, Continued

Methods of Transfection	For established cell lines (e.g., HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. Precisely follow the protocol for your cell line, paying particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in <i>Current Protocols in Molecular Biology</i> (Ausubel <i>et al.</i> , 1994). Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler <i>et al.</i> , 1977), lipid-mediated (Felgner <i>et al.</i> , 1989; Felgner and Ringold, 1989) and electroporation (Chu <i>et al.</i> , 1987; Shigekawa and Dower, 1988). For high efficiency transfection in a broad range of mammalian cells, use Lipofectamine <sup>®</sup> 2000 Reagent available for purchase (see page 12). For more information on Lipofectamine <sup>®</sup> 2000 and other transfection reagents, visit www.lifetechnologies.com or contact Technical Support (see page 13).
Positive Control	pcDNA <sup>TM</sup> 3.1/His/ <i>lacZ</i> is supplied as a positive control vector for mammalian transfection and expression (see page 11). pcDNA <sup>TM</sup> 3.1/His/ <i>lacZ</i> may be used to optimize transfection conditions for your cell line. The gene encoding $\beta$ -galactosidase is expressed in mammalian cells under the CMV promoter. A successful transfection results in $\beta$ -galactosidase expression and can be easily assayed.
Assay for β-galactosidase Activity	You may assay for $\beta$ -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. The $\beta$ -Gal Assay Kit and the $\beta$ -Gal Staining Kit are available for purchase (see page 12 for ordering information) for fast, easy detection of $\beta$ -galactosidase expression.
Geneticin <sup>®</sup> Selective Antibiotic	For stable transfection, pcDNA <sup>™</sup> 3.1/His A, B, and C contain the resistance factor to Geneticin <sup>®</sup> . Geneticin <sup>®</sup> blocks protein synthesis in mammalian cells by interfering with ribosomal function. It is an aminoglycoside, similar in structure to neomycin, gentamycin, and kanamycin. Expression in mammalian cells of the bacterial aminoglycoside phosphotransferase gene (APH), derived from Tn5, results in detoxification of Geneticin <sup>®</sup> Selective Antibiotic (Southern and Berg, 1982).

# Transformation and Transfection, Continued


Geneticin <sup>®</sup> Selection	Geneticin <sup>®</sup> is available for purchase (see page 12 for ordering information). Use as follows:
Guidelines	1. Prepare Geneticin <sup>®</sup> in a buffered solution (e.g., 100 mM HEPES, pH 7.3).
	2. Use 100 to 1,000 $\mu$ g/mL of Geneticin <sup>®</sup> in complete medium.
	3. Calculate concentration based on the amount of active drug (check the lot label).
	<ol> <li>Test varying concentrations of Geneticin<sup>®</sup> on your cell line to determine the concentration that kills your cells (kill curve). Cells differ in their susceptibility to Geneticin<sup>®</sup>.</li> </ol>
	Cells will divide once or twice in the presence of lethal doses of Geneticin <sup>®</sup> Selective Antibiotic, so the effects of the drug take several days to become apparent. Complete selection can take from 2 to 4 weeks of growth in selective medium.
Preparing Cells for Lysis	Use the procedure below to prepare cells for lysis prior to purification of your protein on ProBond <sup><math>M</math></sup> (see page 12 for ordering information). You will need $5 \times 10^6$ to $1 \times 10^7$ cells for purification of your protein on a 2 mL ProBond <sup><math>M</math></sup> column (see ProBond <sup><math>M</math></sup> Protein Purification manual).
	<b>Note:</b> The N-terminal peptide adds approximately 5 kDa to the size of your recombinant fusion protein.
	1. Seed cells in five T-75 flasks or 2 to 3 T-150 flasks.
	2. Grow the cells in selective medium until they are 80–90% confluent.
	3. Harvest the cells by treating with trypsin-EDTA for 2 to 5 minutes or by scraping the cells in PBS.
	<ol> <li>Inactivate the trypsin, if necessary, and transfer the cells to a sterile microcentrifuge tube.</li> </ol>
	5. Centrifuge the cells at approximately $250 \times g$ for 5 minutes. You may lyse the cells immediately or freeze in liquid nitrogen and store at $-80^{\circ}$ C until needed.
Lysing Cells	If you are using ProBond <sup>™</sup> resin, refer to the ProBond <sup>™</sup> Protein Purification manual for details about sample preparation for chromatography. If you are using other metal-chelating resin, refer to the manufacturer's instruction.

#### Appendix

#### Map of pcDNA<sup>™</sup>3.1/His A, B, and C Vectors

Map of pcDNA<sup>™</sup>3.1/His

The figure below summarizes the features of the pcDNA<sup>™</sup>3.1/His vectors. The sequences for pcDNA<sup>™</sup>3.1/His A, B, and C are available for downloading from <u>www.lifetechnologies.com</u> or by contacting Technical Support (see page 13). Details of the multiple cloning sites are shown on page 3 for pcDNA<sup>™</sup>3.1/His A, page 4 for pcDNA<sup>™</sup>3.1/His B, and page 5 for pcDNA<sup>™</sup>3.1/ His C.



# Features of pcDNA<sup>™</sup>3.1/His A, B, and C Vectors

# Features of pcDNA<sup>™</sup>3.1/His

pcDNA<sup>™</sup>3.1/His A (5,514 bp), pcDNA<sup>™</sup>3.1/His B (5,515 bp), and pcDNA<sup>™</sup>3.1/His C (5,513 bp) contain the following elements. All features have been functionally tested.


Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Permits efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987).
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert.
N-terminal polyhistidine tag	Permits purification of your recombinant protein on metal-chelating resin such as ProBond™.
Xpress <sup>™</sup> epitope tag	Allows detection of your recombinant protein with the Anti-Xpress <sup>™</sup> Antibody (Cat. no. R910-25).
Enterokinase cleavage site	Allows removal of the N-terminal polyhistidine tag from your recombinant protein using an enterokinase such as EKMax <sup>™</sup> (Cat. no E180-01).
Multiple cloning site in three reading frames	Allows insertion of your gene and facilitates cloning in frame with the polyhistidine N-terminal tag.
BGH reverse priming site	Permits sequencing through the insert.
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992).
f1 origin	Allows rescue of single-stranded DNA.
SV40 early promoter and origin	Allows efficient, high-level expression of the neomycin resistance gene and episomal replication in cells expressing SV40 large T antigen.
Neomycin resistance gene	Selection of stable transfectants in mammalian cells (Southern and Berg, 1982).
SV40 polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA.
pUC origin	High-copy number replication and growth in <i>E. coli</i> .
Ampicillin resistance gene (β-lactamase)	Selection of vector in <i>E. coli</i> .

## Map of pcDNA<sup>™</sup>3.1/His/*lac*Z

**Description**  $pcDNA^{M}3.1/His/lacZ$  is a 8,577 bp control vector containing the gene for β-galactosidase.  $pcDNA^{M}3.1/His$  B was digested with *Bam*H I, blunted, and digested with *Eco*R I. A 3.2 kb blunt-*Eco*R I fragment containing the β-galactosidase gene was then ligated into  $pcDNA^{M}3.1/His$  B in frame with the N-terminal peptide. The β-galactosidase protein expressed from  $pcDNA^{M}3.1/His/lacZ$  has a size of approximately 120 kDa.

#### Map of Control Vector

The figure below summarizes the features of the pcDNA<sup>™</sup>3.1/His/*lacZ* vector. The nucleotide sequence for pcDNA<sup>™</sup>3.1/His/*lacZ* is available for downloading from <u>www.lifetechnologies.com</u> or by contacting Technical Support (see page 13).



#### **Accessory Products**

#### Additional Products

The following additional products may be used with the pcDNA<sup>™</sup>3.1/His vectors. For more information, visit www.lifetechnologies.com or contact Technical Support (see page 13).

Item	Quantity	Cat. no.
One Shot <sup>®</sup> TOP10 Chemically Competent E. coli	10 reactions	C4040-10
One Shot® TOP10F´ Chemically Competent E. coli	$20 \times 50 \ \mu L$	C3030-03
One Shot <sup>®</sup> Max Efficiency <sup>®</sup> DH5a <sup>™</sup> T1 <sup>R</sup> Competent Cells	$20 \times 50 \ \mu L$	12297-016
MAX Efficiency <sup>®</sup> DH10B <sup>™</sup> Competent Cells	5 × 0.2 mL	18297-010
Electrocomp <sup>TM</sup> Kit	$2 \times 20$ reactions	C66511
Ampicillin	200 mg	11593-027
Carbenicillin	5 g	10177-012
T7 promoter primer	2 µg	N560-02
BGH Reverse primer	2 µg	N575-02
PureLink <sup>®</sup> HiPure Plasmid Miniprep Kit	100 preps	K2100-03
PureLink <sup>®</sup> HiPure Plasmid Midiprep Kit	25 preps	K2100-04
Lipofectamine <sup>®</sup> 2000 Reagent	1.5 mL	11668-019
β–Gal Assay Kit	1 kit	K1455-01
β–Gal Staining Kit	1 kit	K1465-01
Anti-Xpress <sup>™</sup> Antibody	50 µL	R910-25
EKMax <sup>™</sup> Enterokinase	250 U	E180-01
ENVIAX EITTEFORMASE	1000 U	E180-02
Geneticin <sup>®</sup> Selective Antibiotic	1 g	11811-023
	5 g	11811-031
	20 mL (50 mg/mL)	10131-035
	100 mL (50 mg/mL)	10131-027

#### Purifying Fusion Proteins

The N-terminal polyhistidine tag can be used to purify the recombinant fusion protein with a metal-chelating resin such as ProBond<sup>™</sup>. Ordering information for ProBond<sup>™</sup> resin is provided below.

Item	Quantity	Cat. no.
ProBond <sup>™</sup> Purification System	6 purifications	K850-01
ProBond <sup>™</sup> Nickel-Binding Resin	50 mL	R801-01
(Precharged resin provided as a 50% slurry in 20% ethanol)	150 mL	R801-15

# **Technical Support**

Obtaining support	<ul> <li>For the latest services and support information for all locations, go to <u>www.lifetechnologies.com/support</u>.</li> <li>At the website, you can: <ul> <li>Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities</li> <li>Search through frequently asked questions (FAQs)</li> <li>Submit a question directly to Technical Support (<u>techsupport@lifetech.com</u>)</li> <li>Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents</li> <li>Obtain information about customer training</li> <li>Download software updates and patches</li> </ul> </li> </ul>
Safety Data Sheets (SDS)	Safety Data Sheets (SDSs) are available at <u>www.lifetechnologies.com/support</u> .
Certificate of Analysis	The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to <u>www.lifetechnologies.com/support</u> and search for the Certificate of Analysis by product lot number, which is printed on the box.
Limited warranty	Life Technologies Corporation is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about a Life Technologies product or service, contact our Technical Support Representatives. All Life Technologies products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. <u>This warranty</u> <u>limits the Company's liability to only the price of the product</u> . No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order. Life Technologies makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives. Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

#### **Purchaser Notification**

Limited Use Label License No: Research Use Only	The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact <u>outlicensing@lifetech.com</u> or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.
Limited Use Label License 6x His Tag	This product is licensed from Hoffmann-La Roche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

#### References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. J. Biol. Chem. 264, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). Current Protocols in Molecular Biology (New York: Greene Publishing Associates and Wiley-Interscience).
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. Cell 41, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. Molec. Cell. Biol. 7, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. Nucleic Acids Res. *15*, 1311-1326.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. Proc. West. Pharmacol. Soc. 32, 115-121.
- Felgner, P. L. a., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. Nature 337, 387-388.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. J. Biol. Chem. 267, 16330-16334.
- Miller, J. H. (1972). Experiments in Molecular Genetics (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. Molec. Cell. Biol. 7, 4125-4129.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). Molecular Cloning: A Laboratory Manual, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. BioTechniques *6*, 742-751.
- Southern, P. J., and Berg, P. (1982). Transformation of Mammalian Cells to Antibiotic Resistance with a Bacterial Gene Under Control of the SV40 Early Region Promoter. J. Molec. Appl. Gen. 1, 327-339.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. Cell *11*, 223-232.

©2012 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the properties of Life Technologies Corporation or their respective owners.

#### Notes

### Notes

#### Notes



 Headquarters

 5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

 For support visit www.invitrogen.com/support or email techsupport@invitrogen.com



www.lifetechnologies.com